Recent Updates

- ► Logging of target-related parameters
 - ▶ Implemented a general-purpose VI to record a variable array onto TSV
 - ▶▶ In local folder: Desktop/target_data/ and/or
 - ▶▶ **In network folder:** /data2/e1039_data/target_data/
 - Continuous loggings are running fine
 - ▶ Cryo Control
 - ▶▶ Roots-pump LCW
 - ▶▶ Thermocouples for Annealing & Microwave Control Systems
 - Simple data viewer on E1039 data summary
 - https://e906-gat1.fnal.gov/data-summary/e1039/
- Automated alarm via WhatsApp
 - ▶ About NM4 LCW
 - ▶ Will include other systems (after confirming+improving stability)
- Arrangement of GitHub repository
 - Dedicated repositories for Cryo Control and PDP:
 - ▶▶ https://github.com/uva-spin/Cryo-Control
 - ▶▷ https://github.com/uva-spin/PDP
 - ▶ To avoid uncomitted changes and unnecessary conflicts
- ► Measurement of NMR noise level In another presentation